REFRACTOMETRO DIGITAL

Brix Meter

Range 0-32%/55%

Accuracy $\pm 0.2\%$

SOONDA®

REFRACTÓMETRO DIGITAL BRIX

Modelo RD032: Rango de medición 0~32%

Modelo RD055: Rango de medición 0~55%

Características:

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
Measurement steps

01 Power on

02 Drop the sample, at least 2 / 3 of the sample tank

03 Press the "meas" key, results in 5 seconds

04 Clean and dry the sample tank for the next measurement

Cleaning

If there is dirt in the measuring slot or the instrument is not used for a long time, it should be cleaned first.

If the measuring slot contains oiliness or is not easy to be cleaned, please clean it with ethanol or soapy water.

Drip in clean water and wash with paper towels more than three times (throw away the towels after wiping every time)

Zero-setting

Measure whether the clean water or pure water is 0.0% before using, and if it is not 0.0%, it will be zeroed first.

- ① Keep about 1 mm between dripped water or pure water and the top of the measuring slot, and the capacity is about 0.3 ml.

- ② Press "MEAS" key


Press 1 second to zero.

測量
MEAS

- ③ Check if the measured value is 00.0%.

- ④ If the value is not 00.0%, press the zero key to reset.

清零
ZERO

Measurement


Measure again after zero clearing.

- ① Keep about 1 mm between dripped water or pure water and the top of the measuring slot, and the capacity is about 0.3 ml.

- ② Press the "MEAS" key to display the value.

測量
MEAS

Please check the precautions before measuring. Suck out the measuring solution and measure several times to ensure the accuracy and reliability of the data.

Warning

- ▶ This instrument is a precise optical instrument. Prevent it from collision and fall.
- ▶ Forbid hard objects to scratch the optical prism in measuring slot.
- ▶ When it is out of use, please store it in a place as dry as possible.
- ▶ USB charging voltage should not exceed DC6V limit.