


## Vertical Type Shaking Incubator



# About Us

BIOBASE group, the professional manufacturer for Laboratory and Medical equipment. During long term friendly bilateral cooperation relationship with customers, BIOBASE is always being trusted. The customers purchase BIOBASE made goods, at the mean time, more and more customers are willing to buy a variety of goods from BIOBASE.

In order to satisfy customers demand and ensure customers' benefit, BIOBASE made field trip to many suppliers in China, controlled the product quality strictly in stead of customers, as well as strived from the suppliers to get competitive pricing. Recently, BIOBASE has confirmed many OEM partners and will keep searching for more suppliers upon customers' request.

BIOBASE provide one-stop procurement solution which can not only reduce the resourcing cost, but also can ensure good quality and sufficient after-sale service.

In addition, any abroad strong manufacturer is welcomed to contact us, we can work together to open China market through BIOBASE China.



ISO9001



ISO13485



ISO14001



CE



## Large Capacity Vertical Type Shaking Incubator (Double Door & Double Layer)



### Features:

- \* Audio and visual alarm.
- \* Power off recovery function.
- \* Soft shaking mode, stable and reliable operation.
- \* Efficient drive motor and cooling controller, low energy consumption.
- \* Transparent window and built-in lighting design convenient for observation samples.
- \* Accessories: UV lamp(standard accessory); Humidity controller, Inlet air pipe, Light system(optional accessories).



**Technical Parameters:**

Model	BJPX-1102	BJPX-2102
Controller	PID	
Display	LCD	
Circulation Mode	High speed circulating type	
Shaking Mode	Rotary shaking	
Drive Mode	Unishaft balanced drive	
Rotating Speed	30~300rpm	
Rotating Accuracy	±1rpm	
Vibration Range	Φ25mm	
Temp. Range	RT+5~60°C	4~60°C
Temp. Control Accuracy	±0.1°C	
Temp. Uniformity	±1°C(At 37°C)	
Timing Range	0~999.59min	
Shaking Plate Size	734*464 mm	
Shaking Plate Quantity	2 pcs	
Max. Capacity	1000ml*30 or 500ml*44 or 250ml*56 or 50ml*104	
Standard Capacity	500ml*44, 250ml*28	
Internal Material	304 Stainless steel	
External Material	High-quality ductile iron, channel steel, angle steel and cold-rolled steel plate	
Refrigerant	/	CFC-Free, 134A
Cooling Consumption	/	312W
Heating Consumption	600W	
Safety Function	Audio and visual alarm for over-temperature; Alarm for over-speed; Independent over-temperature protector; Shutdown protection for door opening; Leakage or over-current protection; Sensor fault identification; Overload protection of refrigerator	
Consumption	930W	1100W
Power Supply	AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional)	
Internal Size(W*D*H) mm	850*520*760	
External Size(W*D*H) mm	950*740*1400	
Package Size(W*D*H) mm	1140*900*1670	
Gross Weight(kg)	395	



# C E R T I F I C A T E

## ATTESTATION CERTIFICATE OF ELECTROMAGNETIC COMPATIBILITY AND LOW VOLTAGE DIRECTIVES

Technical file of the company mentioned below has been observed  
2014/30/EU Electromagnetic Compatibility Directive and  
2014/35/EU Low Voltage Directives have been taken as references for these processes

Company Name	: Biobase Biodustry (Shandong) Co., Ltd.
Company Address	: No.2717 Gongyeyi Road,Mingshui Economic Developing Zone, Zhangqiu,Jinan,Shandong,China
Related Directives and Annex	: 2014/35/EU Low Voltage Directive 2014/30/EU Electromagnetic Compatibility Directive
Related Standards	: EN 61010-1:2010;EN 61326-1:2013
Product Name	: Incubator
Report No and Date	: OUCE-SH-174058
Product Brand/Model/Type	: BJPX-C50,BJPX-C80,BJPX-C160,BJPX-C270,BJPX-H30,BJPX-H48II,BJPX-H50, BJPX-H64II,BJPX-H80,BJPX-H123II,BJPX-H160,BJPX-H230II,BJPX-H270, BJPX-B80,BJPX-B150,BJPX-B200,BJPX-B250,BJPX-B300,BJPX-B400,BJPX-HT150, BJPX-HT200,BJPX-HT250,BJPX-HT300,BJPX-HT400,BJPX-M150,BJPX-M200, BJPX-M250,BJPX-M300,BJPX-M400,BJPX-A250,BJPX-A300,BJPX-A400, BJPX-A500,BJPX-L150,BJPX-L250,BJPX-L300,BJPX-L400,BJPX-L500,BJPX-G, BJPX-200BX,BJPX-200CX,BJPX-P10,BJPX-P20,BJPX-103B,BJPX-100B,BJPX-200B, BJPX-1102,BJPX-2102,BJPX-1102C,BJPX-2102C,BJPX-2012,BJPX-2012R

Certificate Number	: M.2016.201.N1827
Initial Assessment Date	: 15.08.2017
Registration Date	: 16.08.2017
Reissue Date/No	: -
Expiry Date	: 15.08.2022

  
UDEM International Certification  
Auditing Training Centre Industry  
and Trade Co. Ltd.

The validity of the certificate can be checked through [www.udemltd.com.tr](http://www.udemltd.com.tr). The CE mark shown on the right can only be used under the responsibility of the manufacturer with the completion of EC Declaration of Conformity for all the relevant Directives. This certificate remains the property of UDEM International Certification Auditing Training Centre Industry and Trade Co. Ltd. to whom it must be returned upon request. The above named firm must keep a copy of this certificate for 15 years from the registration of certificate. This certificate only covers the product(s) stated above and UDEM must be noticed in case of any changes on the product(s)  
Address: Mutlukent Mahallesi 2073 Sokak (Eski 93 Sokak) No:10 Çankaya - Ankara - TURKEY  
Phone: +90 0312 443 03 90 Fax: +90 0312 443 03 76  
E-mail: [info@udemltd.com.tr](mailto:info@udemltd.com.tr) [www.udemltd.com.tr](http://www.udemltd.com.tr)



## **BIOBASE GROUP**

Add: No.9 Gangxing Road, High-tech Zone,  
Jinan City, Shandong Province, China  
Tel: +86-531-81219801/03 Fax: +86-531-81219804  
E-mail: [export@biobase.com](mailto:export@biobase.com)  
Web: [www.biobase.cc](http://www.biobase.cc) / [www.meihuatrade.com](http://www.meihuatrade.com)



BIOBASE CHINA



BIOBASE CHINA



BIOBASE CHINA