

JIBI MED[®]

DEVOTE OURSELVES TO THE PROFESSIONAL MANUFACTURING OF THE STERILIZERS

| Company Profile

Jiangyin Binjiang Medical Equipment Co., Ltd.(JIBIMED) is a professional high-tech enterprise, incorporating research, manufacturing, distribution of sterilizing equipment and relevant medical products, located in Jiangyin city of Jiangsu Province, which is a remarkable place with picturesque scenery. Jiangyin is the birthplace of outstanding people, and ranked first on the economy among all the Chinese cities.

With the foundation concept of "technology reflects our wisdom, quality reflects our dignity". JIBIMED has gathered many top professional experts and famous consultants who involving themselves in the research and development of the sterile products after 10 years development. We have formed strong power to follow the tip technology at home and abroad, has become one of the largest and most professional manufacturers with continuous growth of export for the past years.

At present, JIBIMED has two producing zones, covering the area of 40 acres and the building area of 30000 m², we have more than ten full producing line for the vertical pressure steam sterilizer, portable pressure steam sterilizer, table top steam sterilizer, horizontal pressure steam sterilizer, low temperature plasma sterilizer, ultraviolet lamp trolley and intelligent drying oven, reaches 50000 sets for annual production, which are not only covered 70% of the Chinese hospitals, but also delivered to Europe, Africa, East-South Asia, Middle East over 50 countries, we have provided our goods and service to amount of tens of thousands of hospitals all over the world.

All the products from JIBIMED are strictly treated according to quality control system, such as CE, ISO9001, ISO13485. To supply the "zero-defect" units with comprehensive, considerate and sincere service is our final target goal.

The customers satisfaction and the quality expectation is our unswerving pursuit. Advanced technology, scientific management and excellent service is the guarantee of the continuously improving for the quality. Product safety, regulatory compliance and environmental protection is our social responsibility.

JIBIMED has win highly appreciation among domestic and oversea hospitals of various levels in the sterilizing field with innovated products, qualified service, and efficient management, has been awarded the "National Quality Trustworthy Enterprise", "Jiangsu Integrity and Outstanding Enterprises", ect .

For the honor and the commitment , we strive to do better !



Horizontal Pressure Steam Sterilizer



Use

The horizontal autoclave, which adopts the way of gravity exchange to release the cold air from the chamber more completely to ensure the reliable sterilization. The control system automatically adjust the steam inlet and outlet according to chamber temperature during sterilizing. This unit is an ideal equipment for clinics, scientific research institutes and other organizations to sterilize surgical instruments, fabric glassed, and the culture media etc.

Characteristics

- ◆ Automatically controlled the sterilizing process, easy to operate.
- ◆ With drying function, suitable for the medical dressing drying.
- ◆ With over-temperature, over-pressure auto-protect.
- ◆ The door can not be open until the pressure in the chamber is reduced to 0.027MPa. And the unit can not be started on if the door doesn't close well.
- ◆ The safety valve will be open automatically when the inner pressure over 0.24MPa and the steam would be exhausted the water tank.
- ◆ Automatically cut off the power if the water-lack happened, and alarm at the meanwhile.
- ◆ The chamber of the sterilizer is made of stainless steel.

Horizontal Pressure Steam Sterilizer



Technical Data


MODEL	WS-150YDA	WS-200YDA	WS-280YDA	WS-400YDA	WS-500YDA
Technical Data					
Sterilizing chamber volume	150L φ 440 × 1000mm	200L φ 515 × 1000mm	280L φ 600 × 1000mm	400L φ 700 × 1100mm	500L φ 700 × 1300mm
Working pressure	0.22MPa				
Working temperature	134°C				
Temperature adjust range	40~134°C				
Time for sterilization	0~60min				
Time for drying	0~60min				
Heat average	≤ ± 2°C				
Power	9KW/380V.50Hz	9KW/380V.50Hz	12KW/380V.50Hz	18KW/380V.50Hz	18KW/380V.50Hz
Dimension(mm)	1400 × 600 × 1300	1400 × 670 × 1650	1400 × 770 × 1780	1430 × 880 × 1830	1800 × 900 × 1820
Transport dimension(mm)	1560 × 750 × 1850	1560 × 820 × 1850	1680 × 920 × 1950	1600 × 1050 × 2050	1850 × 1050 × 2050
G.W/N.W	320kg / 240kg	350kg / 260kg	465kg / 365kg	530kg / 420kg	580kg / 470kg

High quality of an international norm

As a provincial high-tech enterprise, Binjiang Medical, who relies on the accurate positioning and professional spirit of being brave enough to fight, creates a miracle of development speed in the industry and wins satisfaction and trust from customers.



Professional Accreditation



JIBI MED®

Jiangyin Binjiang Medical Equipment Co., Ltd.

Add: No. 38, Yungu Road, Changshou Community,
Zhouzhuang Town, Jiangyin City, Jiangsu Province

National toll free hotline: **400-700-6180**

Tel: +86-510-86270699 86296768

Fax: +86-510-86270599

P.C.: 214400

Http://www.chengyi-cn.com

E-mail: postmaster@chengyi-cn.com

jibimed@vip.163.com

