



DEVOTE OURSELVES TO THE PROFESSIONAL MANUFACTURING OF THE STERILIZERS

## | Company Profile

Jiangyin Binjiang Medical Equipment Co., Ltd(JIBIMED) is a professional high-tech enterprise, incorporating research, manufacturing, distribution of sterilizing equipment and relevant medical products, located in Jiangyin city of Jiangsu Province, which is a remarkable place with picturesque scenery. Jiangyin is the birthplace of outstanding people, and ranked first on the economy among all the Chinese cities.

With the foundation concept of "technology reflects our wisdom, quality reflects our dignity". JIBIMED has gathered many top professional experts and famous consultants who involving themselves in the research and development of the sterile products after 10 years development. We have formed strong power to follow the tip technology at home and abroad, has become one of the largest and most professional manufacturers with continuous growth of export for the past years.

At present, JIBIMED has two producing zones, covering the area of 40 acres and the building area of 30000 m<sup>2</sup>, we have more than ten full producing line for the vertical pressure steam sterilizer, portable pressure steam sterilizer, table top steam sterilizer, horizontal pressure steam sterilizer, low temperature plasma sterilizer, ultraviolet lamp trolley and intelligent drying oven, reaches 50000 sets for annual production, which are not only covered 70% of the Chinese hospitals, but also delivered to Europe, Africa, East-South Asia, Middle East over 50 countries, we have provided our goods and service to amount of tens of thousands of hospitals all over the world.

All the products from JIBIMED are strictly treated according to quality control system, such as CE, ISO9001, ISO13485. To supply the "zero-defect" units with comprehensive, considerate and sincere service is our final target goal.

The customers satisfaction and the quality expectation is our unswerving pursuit. Advanced technology, scientific management and excellent service is the guarantee of the continuously improving for the quality. Product safety, regulatory compliance and environmental protection is our social responsibility.

JIBIMED has win highly appreciation among domestic and oversea hospitals of various levels in the sterilizing field with innovated products, qualified service, and efficient management, has been awarded the "National Quality Trustworthy Enterprise", "Jiangsu Integrity and Outstanding Enterprises", ect .

For the honor and the commitment , we strive to do better !



## Vertical Pressure Steam Sterilizer



### Use

The vertical pressure steam sterilizers are assembled with heating system, micro-computer controlled system, over heat and over pressure protection system, which are reliable for the sterilizing effects, convenient of operation and energy conservation. They are ideal equipments for clinics, scientific research institutes and other organizations to sterilize surgical instruments, fabrics, glasses, culture media etc.

### Characteristics

- ◆ Fully automatically, no need for supervision, easy to use.
- ◆ Touch key, self-define of the sterilizing data.
- ◆ Hand wheel type of quick-open door structure, self-inflating type seal, and with door safety lock system.
- ◆ Equiped with electric pressure temperature controller, constant voltage and temperature can be automatically kept when the unit reach the setting parameter.
- ◆ Rapid sterilizing, only required 4-6 mins.
- ◆ The sterilizing drum made of high quality stainless steel.
- ◆ Over temperature and over pressure self protection: if the inner pressure is over 0.245MPa, the safety valve will be turn on automatically.
- ◆ The full unit is made of nice quality stainless steel, and stand well.
- ◆ Adopted with immersed heating, with high heating efficiency.
- ◆ Self protection of water lacking.
- ◆ Automatically shut off with beep reminding after sterilization.

### Technical Data

Technical Data	MODEL
	LS-28HD
Chamber volume	28L φ 270 × 500 mm
Working pressure	0.22MPa
Working temperature	134°C
Heat average	≤ ± 1°C
Timer range	0 ~ 99min
Temperature adjust range	105 ~ 134°C
Power	2 KW/AC220V.50HZ
Overall dimension(mm)	440 × 420 × 750
Transport dimension(mm)	490 × 480 × 760
G.W/N.W	28kg / 25kg




# High quality of an international norm

As a provincial high-tech enterprise, Binjiang Medical, who relies on the accurate positioning and professional spirit of being brave enough to fight, creates a miracle of development speed in the industry and wins satisfaction and trust from customers.



# Professional Accreditation





**Jiangyin Binjiang Medical Equipment Co., Ltd.**

Add: No. 38, Yungu Road, Changshou Community,  
Zhouzhuang Town, Jiangyin City, Jiangsu Province

National toll free hotline: **400-700-6180**

Tel: +86-510-86270699 86296768

Fax: +86-510-86270599

P.C.: 214400

Http://www.chengyi-cn.com

E-mail: postmaster@chengyi-cn.com

jibimed@vip.163.com

