

Table Top Pulse Vacuum Steam Sterilizer

TM-XD35/50/90DV



Use

The table type sterilizers are safely auto controlled, which are designed for the dental clinics, hospitals, school lab, scientific research institutions and laboratory. And they are mostly suitable for sterilizing of surgical, dental and eye instruments, glassware medicine, culture medium and biological dressing, food, etc.

Characteristics

- Standard Class B with three times vacuum and drying, the remaining temperature of sterilized instruments less than 0.2%.
- LCD display working status, touch keys.
- The vacuum data can reach -0.8bar, it fits various instrument, including package, unpackaged, solid, hollow, multi apertures and inset pipes. It makes sure the instrument can be sterilized completely.
- Fully computer controlled LCD display bright show. The interface is modularized panel for easy operation.
- It is installed with BOWIE&DICK which measures the penetration of water steam.
- It has vacuum test which can make sure the ability of vacuum testing. Fast and independent steam generator.

- Built-in mini printer which can record the sterilizing information.
- Equipped with Safety valve, Safety lock door system, Pressure lock system, Pressure or temperature overload protection, Alarm system.
- It inactivates the HIV, HBV, BSE, and bacillus efficiently.

Technical Data

MODEL	TM-XD35DV	TM-XD50DV	TM-XD90DV
Technical Data			
Sterilizing chamber volume	35L φ300×500mm	50L φ340×550mm	90L φ400×720mm
Working pressure	0.22Mpa		
Working temperature	134℃		
Temperature adjust range	105-134℃		
Timer range	0-99min		
Heat average	≤±1℃		
Power	AC110/220V 60Hz/2.7KW		AC110/220V 60Hz/4KW
Dimension (mm)	730×550×510	770×600×540	1000×690×610
Exterior package dimension (mm)	810×680×650	880×730×700	1070×740×700
G.W/N.W	115kg/93Kg	130kg/106Kg	200kg/180Kg

Details

