



BIOBASE®

Version 2021-1

About Us

Founded in 1999, having over 2000 employees (till August, 2018), BIOBASE Group is a new high-tech enterprise focusing on development, production and management of scientific equipment, In Vitro Diagnosis (IVD) instrument and reagent.

BIOBASE Group is specialized in products of 8 areas including medical diagnosis, biosafety protection, disinfection and sterilization solution, water purification system, infant care products, cold chain products, software products, clean room project.

We have full production line for table top steam sterilizer, vertical pressure steam sterilizer, portable pressure steam sterilizer, horizontal pressure steam sterilizer, hot air sterilizer, gas sterilizer, glass bead sterilizer, infrared sterilizer and automatic glassware washer. And every year we have many new products that come into the market.

BIOBASE Group has certificates including ISO9001, ISO13485, ISO14001, SGS, CE, etc for various products.

BIOBASE has set up branch offices in 24 provinces in China, and has also set up warehouses in the United States, India, UAE and Jordan. We now have long-term business partners in more than 130 countries.

BIOBASE Group always adheres to the mission of "service for health for the future" and promotes company's culture vigorously. Providing with high quality products and good service, BIOBASE would like to establish win-win business cooperation with distributors all over the world.



ISO9001



ISO13485



ISO14001



CE



Toxic Chemical Storage Cabinet



CSC-12D



CSC-22D



CSC-30D



CSC-45D



CSC-60D



CSC-90D

Features:

- * Two layers of fireproof steel plates with 38mm thermal insulation between the two steel plates.
- * Main body material: 1.0mm steel plate, epoxy resin spraying.
- * Security: Equipped with combination lock and mechanical lock.
- * 50mm high leak-proof liquid tank ensures no liquid spillage.
- * Installed with dual air outlets with fire protection.
- * Customizable.

Technical Parameters:

Model	CSC-12D	CSC-22D	CSC-30D	CSC-45D	CSC-60D	CSC-90D
Capacity	12gallon (45L)	22gallon (83L)	30gallon (114L)	45gallon (170L)	60gallon (227L)	90gallon (340L)
Standard Shelf	1pc	3pcs	1pc	2pcs	2pcs	2pcs
External Size (W*D*H)mm	590*460*890	600*460*1650	1090*460*1120	1090*460*1650	860*860*1650	1090*860*1650
Package Size (W*D*H)mm	665*535*915	680*540*1795	1170*540*1260	1170*540*1790	1790*940*940	1790*940*1790
Gross Weight	50kg	93kg	101kg	138kg	170kg	200kg

BIOBASE®

Add: No.9 Gangxing Road, High-tech Zone, Jinan City,
shandong Province, China

Tel: +86-531-81219803/01 Fax: +86-531-81219804

Email: export@biobase.com

Website: www.biobase.cc www.meihuatrade.com



BIOBASE CHINA



BIOBASE CHINA



BIOBASE CHINA