

**ACCUMAX**

**FINESSE &  
CLARITY  
PERSONIFIED**

**ACCUMAX**

**PETRI**

**PLATES**

## ACCUMAX PETRI PLATES

Accumax produces high-quality Petri dish which is manufactured under strict aseptic conditions with USP class VI complaint polystyrene.

Designed specifically for microbiology applications. With a 90 mm diameter, they are perfect for use in automated plate filing systems. Vented and non-vented variants are available.

Vented dishes feature an elevated lid, promoting excellent gas exchange and evaporation. This design allows for optimal conditions for growth and incubation. On the other hand, non-vented dishes have a flat lid that prevents evaporation, encouraging a longer incubation period and facilitating anaerobic applications.



## Product Features



Manufactured with USFDA  
21 CFR compliant polystyrene



One compartment  
with both vented and  
non vented options



Fully automated  
manufacturing and packaging



Endotoxin/Pyrogen free

Petri plates with  
90 mm diameter



Produced in  
Class 10K clean rooms



Raw materials are compliant  
with ISO 1133 & 306 guidelines



Sterility Assurance level of  
Sterile plates is  $10^{-6}$  and that  
of Aseptic plates is  $10^{-3}$  SAL



## Product Features



Fully stackable with raised ridge and optimally clear petri plate

Petri plates with excellent mechanical and heat resistance properties



Two Packaging options: 10 plates/ bag & 20 plates/ bag

Lot number, Sterility & expiration traceability in bag packaging



These Raised Ridge plates are designed for easy stacking by hand. They stack high and stay stacked without slipping, sliding, or scratching.

These plates are fully compatible with medium throughput carousel based media preparation systems.



# Ordering Information

## Petri Plates

Cat No.	Product Description	Packing Type	Quantity / bag	Quantity bag / Case	Total Quantity per Case
AP-90V-10	Petriplates 90 mm diameter, Vented, Aseptic, Raised Ridge	Bulk	10	50	500
AP-90NV-10	Petriplates 90 mm diameter, Non Vented, Aseptic, Raised Ridge	Bulk	10	50	500
AP-90V-20	Petriplates 90 mm diameter, Vented, Aseptic, Raised Ridge	Bulk	20	25	500
AP-90NV-20	Petriplates 90 mm diameter, Non Vented, Aseptic, Raised Ridge	Bulk	20	25	500
AP-90V-10-S	Petriplates 90 mm diameter, Stack Size 10, Vented, Sterile, Raised Ridge	Bulk	10	50	500
AP-90NV-10-S	Petriplates 90 mm diameter, Stack Size 10, Non Vented, Sterile, Raised Ridge	Bulk	10	50	500
AP-90V-20-S	Petriplates 90 mm diameter, Stack Size 20, Vented, Sterile, Raised Ridge	Bulk	20	25	500
AP-90NV-20-S	Petriplates 90 mm diameter, Stack Size 20, Non Vented, Sterile, Raised Ridge	Bulk	20	25	500

# ACCUMAX

Accumax Lab Devices Pvt. Ltd.  
www.accumaximum.com



ISO 9001:2015  
ISO 13485:2016

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.