



**BIOBASE®**

Version 2021-1

# About Us

Founded in 1999, having over 2000 employees (till August, 2018), BIOBASE Group is a new high-tech enterprise focusing on development, production and management of scientific equipment, In Vitro Diagnosis (IVD) instrument and reagent.

BIOBASE Group is specialized in products of 8 areas including medical diagnosis, biosafety protection, disinfection and sterilization solution, water purification system, infant care products, cold chain products, software products, clean room project.

We have full production line for table top steam sterilizer, vertical pressure steam sterilizer, portable pressure steam sterilizer, horizontal pressure steam sterilizer, hot air sterilizer, gas sterilizer, glass bead sterilizer, infrared sterilizer and automatic glassware washer. And every year we have many new products that come into the market.

BIOBASE Group has certificates including ISO9001, ISO13485, ISO14001, SGS, CE, etc for various products.

BIOBASE has set up branch offices in 24 provinces in China, and has also set up warehouses in the United States, India, UAE and Jordan. We now have long-term business partners in more than 130 countries.

BIOBASE Group always adheres to the mission of "service for health for the future" and promotes company's culture vigorously. Providing with high quality products and good service, BIOBASE would like to establish win-win business cooperation with distributors all over the world.



ISO9001



ISO13485



ISO14001



CE



## Gas Cylinder Cabinet



BKSC-1Q



BKSC-2Q



BKSC-3Q

### Introduction:

The gas cylinder cabinet is a metal container used to improve local exhaust ventilation, protect steel cylinders (gas cylinders) from fire outside the cabinet and protect surrounding objects from internal fire.

### Features:

- \* Cabinet body: It adopts high-quality cold-rolled steel plate, which is treated with pickling and phosphating, and the surface is electrostatically sprayed to achieve the effect of anti-acid and alkali and anti-rust.
- \* Door panel: It adopts detachable hinges, and the window is tempered glass.
- \* Ventilation holes: There are 4 ventilation holes on each side to facilitate gas circulation.
- \* Fixed chain: Internal fixed chain is used to prevent the cylinder from tipping over.
- \* Pedal: There are adjustable pedals at the bottom of the cabinet to facilitate the loading and unloading of gas cylinders.

### Technical Parameters:

Model	BKSC-1Q	BKSC-2Q	BKSC-3Q
Capacity	40L(1 bottle)	80L(2 bottles)	120L(3 bottles)
Door Type	Single door	Double door	
External Size(W*D*H)	600*450*1900mm	900*450*1900mm	1200*450*1900mm
Package Size(W*D*H)	674*524*2030mm	974*524*2030mm	1274*524*2030mm
Net Weight	61kg	85kg	102kg
Gross Weight	78kg	105kg	124kg

### Optional Alarm Apparatus:

NO.	NO.1	NO.2
Power Supply	220V 50Hz	110~220V 50/60Hz
Alarm Type	Audible and visible alarm	
Identifiable Gas	methane, acetylene, hydrogen, coal gas	methane, butane, ethanol, hydrogen, liquefied gas, acetylene, ammonia, etc.

# BIOBASE<sup>®</sup>

Add: No.9 Gangxing Road, High-tech Zone, Jinan City,  
shandong Province, China

Tel: +86-531-81219803/01 Fax: +86-531-81219804

Email: [export@biobase.com](mailto:export@biobase.com)

Website: [www.biobase.cc](http://www.biobase.cc) [www.meihuatrade.com](http://www.meihuatrade.com)



BIOBASE CHINA



BIOBASE CHINA



BIOBASE CHINA