



SKE801

NURSE STOOL



SAIKANG / Medical

Healthier
your life

SKE801

FEATURES AND HIGHLIGHTS

Technical parameters

- Seat diameter ϕ 360mm
- Inclusive wheel height 450-570mm



Structural Design

Freely adjustable up and down, meeting the height needs of different people, 360° rotation to the left and right for different work scenarios.



Stool Surface

Covered with artificial leather, with native cotton, density 45-55D.



Gas Spring

Height adjustment with gas spring, controlled by handle.




Chassis

The base is chrome-plated, designed ergonomically, allowing feet to be more comfortably flat, sitting more comfortably.




Castors

Made of PA nylon material, easy to push and light, steering smooth and wear-resistant.





CONSTANT INNOVATION




Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical (Philippines) Company Limited

Manila, The Philippines
TEL: +63 953 044 7864

US Branch

Healthward Technology, INC

17800 CASTLETON ST STE 180,
CITY OF INDUSTRY CA 91748
TEL: +1 (603)416-3063

Russia Branch

LLC 'SAIKANG'

454004, Russian Federation,
Chelyabinsk,
Universitetskaya
Naberezhnaya, Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC

Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ export@saikangmedical.com

🌐 www.saikangmedical.com



Company video



General catalogue