



# X07a

MESA DE EXAMEN

SK-A1



**SAIKANG** / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE  
Stock Code: 870098

## Parámetros técnicos

L\*W\*H

1900\*600\*680mm

## Parámetros de configuración

Cubrepiés Anti-slip

4pcs

## Piezas ricas



1.2-1.5mm Espesor del bastidor de acero con recubrimiento de polvo 250kg de carga soportada



4 fundas antideslizantes para pies de cama



## Control de calidad

### Cuero

Resistente al fuego, resistente al moho, resistente a la flexión, 100.000 veces resistente a la intemperie

100,000<sub>times</sub>



### Esponja

4cm de espesor, espuma de alta densidad. Densidad 30kg/m

4<sub>cm</sub> | 30<sub>kg/m</sub>

### Pintura

Proceso de pintura epoxi. ASTM pruebas anti-material, espesor de la pintura 0,12mm. Brillo de 60, la pintura puede resistente de 50kg de impacto

0.12<sub>mm</sub> | 60° | 50<sub>kg</sub>



### Material metalico

Baosteel en la lista fortune Global 500 Panasonic Robotic garantiza un acuñado suave a 360

500 | 360°







## INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

### JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901   
[www.chinasaikang.com](http://www.chinasaikang.com)

#### SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,  
Jiangsu Province,China.  
TEL: +86-512-58525850  
FAX: +86-512-58520713  
E-mail: [inquiry@jssaikang.com](mailto:inquiry@jssaikang.com)

#### SAIKANG (USA Branch)

Saikang Medical Technology, Inc  
17800 Castleton St Ste 180,  
City of Industry, Ca 91748  
Tel: +1 949648-2776

#### SAIKANG (UK Branch)

Uk Healthward International Limited  
Owner(s) Address DEPT 114,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +44 2032399738

#### SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.  
Viraj Complex Unit 3 Mombasa Road,  
PO Box 14262-00800, Nairobi.  
TEL: +254 733-704 690

#### SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC.  
X22-S03-England Cluster,  
International City, Dubai UAE  
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.  
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.