


BALANCES



VERSION.2022

ABOUT US

BIOVOPEAK is a research-and service-driven enterprise in laboratory field with the mission of continuously improving the intelligence, precision, safety and convenience of the laboratory. Based on independent manufacturing capability, professional integration of worldwide laboratory resources and localized after-sales service network, we provide one-stop service for all customers.



0.01g Precision Balance

BP-112-2C BP-212-2C BP-312-2C BP-512-2C
BP-612-2C BP-1102-2C BP-2102-2C BP-3102-2C



Specifications

Model	BP-112-2C	BP-212-2C	BP-312-2C	BP-512-2C
Max. capacity(g)	110	210	310g	510
Cal.weight	100	200	200	200
Readability(g)	0.01			
Repeat ability(g)	0.01			
Linearity(g)	0.01			
Stable Time	≤3s			
Pan Size	Φ130mm、160*160mm			
Operation temp	5-25℃			
Option	Printer,Draft shield,RS232,Weight			
Outline size	270x195x85mm			
Unit set Packing Size	310x230x130mm			
Ten Packing Size	660x470x320mm			
N.weight	Unit set 2.5kg/Ten 25kg			
G.weight	Unit set 3kg/Ten30kg			

Description



LCD/AC and DC power supply



Tare function/Counting/Unit conversion (g/ct/oz)...



Min weighing set



Overload alarm/Level indicator



Option: Dual display/Interface/Printer/Windproof/Dust cover

BP-612-2C	BP-1102-2C	BP-2102-2C	BP-3102-2C
610	1100	2100	3100
200	500	Optional	Optional
0.01			
0.01			
0.01			
≤3s			
Φ130mm、160*160mm			
5-25℃			
Printer,Draft shield,RS232,Weight			
270x195x85mm			
310x230x130mm			
660x470x320mm			
Unit set 2.5kg/Ten 25kg			
Unit set 3kg/Ten30kg			



LOW VOLTAGE DIRECTIVE, ELECTROMAGNETIC COMPATIBILITY DIRECTIVE ATTESTATION OF CONFORMITY

Technical file of the company mentioned below has been inspected and audit has been completed successfully.

2014/30/EU Electromagnetic Compatibility Directive and 2014/ 35/EU Low Voltage Directive has been taken as referances for these processes.

Company Name : **Infitek Co., Ltd.**

Company Address : Rm. 2014, Bldg. 3, LigaoguojiHuayuan, No. 1222, West Aoti Road, Lixia District, Jinan, Shandong, China

Related Directives and Annex : **2014/35/EU Low Voltage Directive /Annex III**
2014/30/EU Electromagnetic Compatibility Directive /Annex II

Related Standards : **EN 61010-1:2010/A1:2019, EN IEC 61326-1:2021**

Product Name : **Balance**

Report No and Date : EC.INFI.20230625-R21

Product Brand/Model/Type : BMA Series, BP series, BSP series, BS-WP series, BA Series, BA-B series, BA-C series, BA-D series, BA-E series, BA-F series, BA-G series, BA-H series, BA-M series, BA-N series, BA-P series, BA-R series, BA-S series, BA-T series, BA-V series, BA-W series, BA-Y series, HB Series, MCA Series

Certificate Number : **M.2023.206.C87062**

Initial Assessment Date : 29.06.2023

Registration Date : 30.06.2023

Reissue Date/No : -

Expiry Date : **29.06.2028**

UDEM International Certification
Auditing Training Centre Industry
and Trade Inc. Co.

The validity of the certificate can be checked through www.udem.com.tr. Upon completion of EC declaration of conformity, it is used solely at the manufacturer's responsibility. This certificate remains the property of UDEM International Certification Auditing Training Centre Industry and Trade Inc. Co. to whom it must be returned upon request. The above named firm must keep a copy of this certificate for 15 years from the registration of certificate. This certificate only covers the product(s) stated above and UDEM must be noticed in case of any changes on the product(s)

Address: Mutlukent Mahallesi 2073 Sokak (Eski 93 Sokak) No:10 Çankaya – Ankara – TURKEY

Phone: +90 0312 443 03 90 **Fax:** +90 0312 443 03 76

E-mail: info@udem.com.tr www.udem.com.tr



FDA

Fiscal Year 2023

CERTIFICATION OF REGISTRATION

This certifies that:

Name: BIOEVOPEAK CO., LTD.

Add: Rm. 2014, Bldg. 3, Ligaoguojiuayuan, No. 1222, West Aoti Road, Lixia District, Jinan, Shandong, 250100, CHINA

has completed the FDA Establishment Registration (as manufacturer and foreign exporter) and Device Listing with the US Food & Drug Administration, through Registration/FEI Number: 3020526004

Owner/Operator Number: 10083884

ABmed will confirm that such registration remains effective upon request and presentation of this certificate until the end of the year stated above, unless said registration is terminated after issuance of this certificate. ABmed makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate - holder's device or establishment by the U.S Food and Drug Administration. ABmed assumes no liability to any person or entity in connection with foregoing.

Date of verification: Dec.29, 2022

Date of expiration: Dec.31, 2023

SH OFFICE

TEL:0086-21-50313932 Boyle Wang Phone:0086-18930777676 info@truthful.com.cn

ABMED SERVICE INC.

1312 17th Street Suite 692 Denver, CO 80202

TEL:303-800-0162 info@abmed.com.cn

FDA

Listing No	Premarket Submission NO.	Code	Device Name
D462801	510(k) exempt	DOP	MASS SPECTROMETER, CLINICAL USE, COD Analyzer, Conductivity Meters, Dissolved Oxygen Meters, Ion Meters, Multiparameter Meters, ORP Meters, Benchtop PH Meter, Pocket PH Meter, Salinity Meters, TDS Meters, Titrator, Water Hardness Meters, Chlorine Meter
D462814	510(k) exempt	JJN	APPARATUS, ELECTROPHORESIS, FOR CLINICAL USE, Electrophoresis Cell, Nucleic Acid Sequencing, Electrophoresis Cell, Vertical Type, Electrophoresis Tank, 2-dimensional, Electrophoresis Tank, Isoelectric Focusing Type, Electrophoresis Tank, Semi-dry Transfer Type, Horizontal Electrophoresis Cell, Inner Core of Electrophoresis Tank, Sequencing Electrophoresis Tank, Electrophoresis Power Supply, Electrophoresis System, Gel Imaging System, Transilluminator,
D462815	510(k) exempt	JJQ	COLORIMETER, PHOTOMETER, SPECTROPHOTOMETER FOR CLINICAL USE, Turbidity Meters, Colorimeter, Gloss Meter, Hazemeter, Melting Point Apparatus, Polarimeters, ABBE Refractometer, Automatic Refractometer, Portable Refractometer, Atomic Absorption Spectrophotometer, Flame Photometer, Flame Photometer, Grating Spectrophotometer, Micro-Volume UV/VIS Spectrophotometer, NIR Spectrophotometer, Scanning Spectrophotometer, Single Beam Spectrophotometer, UV Visible Spectrophotometer
D462811	510(k) exempt	IDL	MICROTOME, ACCESSORIES, Digital Pathology Slide Scanner, Paraffin Dispenser, Paraffin Trimmer, Slide Cabinet, Tissue Flotation Water Bath, Tissue Processor

FDA

Listing No	Premarket Submission NO.	Code	Device Name
D462826	510(k) exempt	LXG	EQUIPMENT, LABORATORY, GENERAL PURPOSE, LABELED OR PROMOTED FOR A SPECIFIC MEDICAL USE, Isolator, IVC Cage, IVC System, Analytical Balances, Moisture Balances, Platform Scales, Precision Balances, Semi-Micro Balance, Acid and Alkali Storage Cabinet, Chemical Safety Cabinet, Gas Cylinder Storage Cabinet, Hospital Furniture, Bullet Ice Maker, automatic type, Snow Ice Maker, automatic, CO2 Incubators, CO2 Stacked Thermostatic Oscillator, Constant Temp. & Humidity Chamber, Cooling Biochemical Incubator, BOD Incubator, Cooling Biochemical Incubator, Microbiology Cooling Incubator, Egg Incubator, Heating Incubator, Plant Growth Chamber, Shaking Incubators, Glassware Washers, Ultrasonic Cleaners, Biological Safety Cabinets, Dispensing Booth, Fume Hood, Laminar Flow Cabinet, Liquid Dispenser, Liquid Vacuum Aspiration, Pipettes, Density Meter, Viscosity Baths, Digital Viscometer, Total Organic Carbon Analyzer, Microplate Reader, Microplate Washer, Block Heater, In Situ Hybridization, Thermal Cycler, Clarification Tester, Transparency Tester, Disintegration Testers, Gelatin Testing, Tablet Friability Testers, Tablet Hardness Testers, Tablet Thickness Testers, Medicine Stability Testing Chamber, Homogenizers, Fully Automatic Sample Digestion System, Graphite Digester, Microwave Digestion System, Microwave Digestion/Extraction System, Microwave Sample Preparation System, Sample Concentrators,

FDA

Listing No	Premarket Submission NO.	Code	Device Name
D462823	510(k) exempt	KSQ	DEVICE, BLOOD MIXING AND BLOOD WEIGHING,Blood Collection Monitor,Blood Plasma Extractor
D462427	510(k) exempt	KZJ	DEVICE, GAS GENERATING,Anaerobic Incubators,Anaerobic Jar,
D462824	510(k) exempt	KZA	DEVICE, VEIN LOCATION, LIQUID CRYSTAL,Vein Finder,
D462821	510(k) exempt	JTQ	BATH, INCUBATORS/WATER, ALL,Thermostatic Magnetic Stirring Bath.Water Bath Oscillator, Thermostatic,Water Bath Shaker, TFT Display,Water Bath,Water Bath,Water bath, WB8T
D462822	510(k) exempt	KSE	REFRIGERATOR, FREEZER, BLOOD STORAGE,Freezers, Refrigerator,Blood Bank Refrigerator,Car Refrigerator,Combined Refrigerator & Freezer,Corpse Refrigerator,Portable Refrigerator
D462812	510(k) exempt	IDP	MICROTOME, CRYOSTAT, Floation Workstation, Cryostat Microtome, Fast Freezing Machine, Fully Automated Microtome, Manual Rotary Microtome, Rotary Microtome, Semi-Automated Rotary Microtome, Vibrating Microtome
D462813	510(k) exempt	JDY	EVACUATOR, VAPOR, CEMENT MONOMER,Filter Funnel ManifoldsChemistry-design Diaphragm Pump,Diaphragm Pump,Diaphragm Vacuum Pump,Water Circulating Vacuum Pump
D462810	510(k) exempt	IBJ	LIGHT, MICROSCOPE,Biological Microscope,Digital Microscope,Fluorescence Microscope,Inverted Microscope,Metallurgical Microscope,Polarizing Microscope,Stereoscopic Microscope,Teaching Microscope

FDA

Listing No	Premarket Submission NO.	Code	Device Name
D462820	510(k) exempt	JRS	PURIFIER, WATER (ABSORPTION, DEIONIZATION, MEMBRANE FILTER, REVERSE OSMOSIS),Water Distiller,Double-stage RO Pure Water Machine,Ultra-pure Water Machine
D462818	510(k) exempt	JRQ	SHAKER/STIRRER,Combination Stirring Hotplates,Conventional Hotplates,Heating Mantles,Magnetic Stirrers,Overhead Stirrer,Roller Mixer,Rotary Mixer,Vortex Mixer,Linear Shaker,Orbital Shaker,Rocker Shaker,
D462819	510(k) exempt	JRR	REGULATOR, TEMPERATURE,Cooling Circulator,Low Temperature Circulator,Recirculating Chiller,Refrigerated and Heating Circulator,20L Industrial Digital Rotary Evaporator,Dynamic Temperature Control System,Vacuum Controller,
D462816	510(k) exempt	JQC	CENTRIFUGES (MICRO, ULTRA, REFRIGERATED) FOR CLINICAL USE,Clinical Centrifuge, PRP,Fat Transplantation Centrifuge,Gel Card Centrifuge,Microplate Centrifuge, PCR,Oil Test Centrifuge,PRP Centrifuge,High Speed Centrifuges,High Speed Refrigerated Centrifuges,Low Speed Centrifuges,Low Speed Refrigerated Centrifuges,Mini Centrifuges,Multi Purpose Centrifuges
D462817	510(k) exempt	JRJ	UNIT, DRYING,Dual Convection Oven,Dual Purpose (Oven/Incubator),Forced-Air Drying Oven,Natural Convection Oven,Vacuum Oven,
D462809	510(k) exempt	FLH	WASHER OF BODY WASTE RECEPTACLES, Eye Washer
D462825	510(k) exempt	LDM	INSTRUMENTATION, HIGH PRESSURE LIQUID CHROMATOGRAPHY, Liquid Chromatograph

FDA

Listing No	Premarket Submission NO.	Code	Device Name
D474045	510(k) exempt	GKB	DEVICE, AUTOMATED SEDIMENTATION RATE - ESR Analyzer
D474049	Enforcement Discretion	CHL	Electrode measurement, blood-gases (pco2, po2) and blood ph - Blood Gas Analyzer
D474046	510(k) exempt	KQO	AUTOMATED URINALYSIS SYSTEM - Urine Analyzer
D474047	510(k) exempt	GGN	Plasma, coagulation control - Coagulation Analyzer
D474048	510(k) exempt	JJP	ELECTRODE, ION SELECTIVE (NON-SPECIFIED) - Electrolyte Analyzer
D474050	Enforcement Discretion	FAZ	SYSTEM, CRYOSURGICAL, LIQUID NITROGEN, FOR UROLOGY - Liquid Nitrogen Tank



Bioevopeak

TEL: +86-531-88982330

FAX: +86-531-88983691

Website: Bioevopeak.com

Email: info@bioevopeak.com

Service: support@bioevopeak.com

Address: 17th Floor, Mingsheng Building, High-tech Zone, Jinan City