



BIOBASE®



BIOBASE[®]

INTRODUCTION

BIOBASE Group is a professional manufacturer of laboratory products since 1999. Our main products are Biosafety Cabinet (Class I / II / III), Laminar Flow Cabinet, Fume Hood, Biochemistry Analyzer, Reagent, Autoclave and Incubator. We have more than 2000 employees and 20000 square meters production base. Our output per year reaches 10000 units. Pursuing the highest quality is the tenet of BIOBASE. We hope everybody will enjoy BIOBASE service to keep away from the infectious disease.

Our advantages:

- We built an excellent research team which consists of 50 professional engineers.
- We have 20 years production experience and a strong technical production team. And we adopt high-tech and state-of-the-art production equipment from Italy, Japan, and Switzerland, etc.
- We have a strong quality control team to inspect all production stage, from the purchased accessories to the finished products.
- We have full sets of advanced quality assurance equipments.
- We have ISO9001 quality management system; ISO14001 environment management system, ISO3485, CE mark, MET and NSF 49 certificate.
- We got EN 12469 certificate for Biosafety Cabinet in 2013.




Fume Hood


Fume Hood is used to protect lab environment and operator during general chemical applications. It actively protects operator from inhaling toxic vapors and dramatically reduces the risk of fire and explosion. By installing proper filter, it can also protect environment.

**Ductless
Fume Hood**






Ductless PP Fume Hood



FH700(PD)



FH1000(PD)

Advantage:

- * The case is made of 8mmPP plate, with good finish and corrosion resistance.
- * The front window has 8°physiological incline, ergonomic design, larger viewing angle.
- * The front window is foldable acrylic window, and the lower acrylic window has the function of stopping at will.
- * Electronic control system with anti-overload, anti-shock and other functions, stable performance, long service life.
- * With power failure memory function, that is, when encountering a sudden power failure, power again can maintain the working state before the power failure, convenient experimental operation.

Technical Parameters:

Model	FH700(PD)	FH1000(PD)
External Size	700*700*1100mm	1000*700*1100mm
Internal Size	680*606*688mm	980*606*688mm
Max Opening	640mm	
Airflow Volume	230 m ³ /h	
Air Velocity	0.3~0.8m/s	
Noise	≤65dB(A)	
LED Lamp	15W	
Illumination	≥400lx	
Blower	Built-in centrifugal blower	
Power Supply	220V, 50/60Hz	
Consumption	150W	
Material	Main Body: PP plate, thickness 8mm Work Table: Chemical resistant phenolic resin	
Standard Accessory	Active carbon filter*2	
Optional Accessory	Water tap, Gas tap, Water sink, Base stand, UV lamp, HEPA filter	
Net Weight	67kg	85kg
Gross Weight	87kg	115kg
Package Size	820*820*1280mm	1120*820*1280mm



REGISTRATION NO. 04722Q10272R2M

CERTIFICATE OF QUALITY MANAGEMENT SYSTEM

This is to certify that the quality management system of

Jinan Biobase Biotech Co., Ltd.

Registered Address: OLABO Intelligent Manufacturing Industrial Park, Ancheng Town, Pingyin County, Jinan City, Shandong

Manufacturing Address: North Side of Jiwang Road, Mingshui Economic Development District, Zhangqiu, Jinan City; OLABO Intelligent Manufacturing Industrial Park, Ancheng Town, Pingyin County, Jinan City, Shandong; Workshop No.7, Medical Device Industrial Park, Juancheng County, Heze City, Shandong Province

**Has been assessed and conformed to the following standard(s)
GB/T 19001-2016 idt ISO 9001:2015**

The certificate is valid for the following scope:

The design, development, production and service of Biological Safety Cabinet; Class I Biosafety Cabinet; CO₂ Incubator; Low Temperature Freezer; Medical Refrigerator; Laminar Flow Cabinet; Biological Isolation Chamber; Fume Hood; PCR Cabinet; Dispensing Booth.

Date of issue: July 22, 2022

Date of expiry: July 21, 2025

Date of change: October 20, 2022

General Manager:

**BEIJING HUA GUANG CERTIFICATION
OF MEDICAL DEVICES CO., LTD.**



中国认可
国际互认
管理体系
MANAGEMENT SYSTEM
CNAS C047-M

Note: This certificate will not be continuously valid until the organization has been approved in the annual surveillance audit. The certificate information are available on the website of the certification and accreditation administration of the People's Republic of China (<http://www.cnca.gov.cn>) or the website of CMD (<http://www.cmdc.com.cn>). Address: 5th floor of Zhong Lian building, No. jia88, An Ding Men Wai street, Dongcheng district, Beijing, 100011, P.R. China Telephone: 010-62351993



Registration No. 04722E10003R0M

CERTIFICATE OF ENVIRONMENTAL MANAGEMENT SYSTEM

This is to certify that:

Jinan Biobase Biotech Co., Ltd.

Registered address: **OLABO Intelligent Manufacturing Industrial Park, Ancheng Town, Pingyin County, Jinan City, Shandong**

Manufacturing address: **North Side of Jiwang Road, Mingshui Economic Development District, Zhangqiu, Jinan City; OLABO Intelligent Manufacturing Industrial Park, Ancheng Town, Pingyin County, Jinan City, Shandong**

Operates an Environmental Management System which complies with the requirements of:

GB/T 24001-2016 idt ISO 14001:2015

For the following scope:

Environmental management activities of design, development, manufacture, distribution and service of Biological Safety Cabinet, Laminar Flow Cabinet, Biological Isolation Chamber, Air Purifier, Fume Hood, PCR Cabinet, Dispensing Booth.


Date of issue: **June 14, 2022**

Date of expiry: **June 13, 2025**

Date of change: **June 04, 2023**

General manager:

**BEIJING HUA GUANG CERTIFICATION
OF MEDICAL DEVICES CO., LTD.**



中国认可
国际互认
管理体系
MANAGEMENT SYSTEM
CNAS C047-M

Note: This certificate will not be continuously valid until the organization has been approved in the annual surveillance audit. The certificate information is available on the website of the certification and accreditation administration of the People's Republic of China (<http://www.cnca.gov.cn>) or the website of CMD (<http://www.cmdc.com.cn>). Address: 5th floor of Zhong Lian building, No. jia88, An Ding Men Wai street, Dongcheng district, Beijing, 100011, P.R. China Telephone: 010-62351993



REGISTRATION NO. 04722Q10000272

CERTIFICATE OF QUALITY MANAGEMENT SYSTEM FOR MEDICAL DEVICES

This is to certify that the quality management system of

Jinan Biobase Biotech Co., Ltd.

Registered Address: OLABO Intelligent Manufacturing Industrial Park, Ancheng Town, Pingyin County, Jinan City, Shandong

Manufacturing Address: North Side of Jiwang Road, Mingshui Economic Development District, Zhangqiu, Jinan City; OLABO Intelligent Manufacturing Industrial Park, Ancheng Town, Pingyin County, Jinan City, Shandong; Workshop No.7, Medical Device Industrial Park, Juancheng County, Heze City, Shandong Province

**Has been assessed and conformed to the following standard(s)
YY/T 0287-2017 idt ISO 13485:2016**

The certificate is valid for the following scope:

The design, development, production and service of Biological Safety Cabinet; Class I Biosafety Cabinet; CO₂ Incubator; Low Temperature Freezer; Medical Refrigerator; Laminar Flow Cabinet; Biological Isolation Chamber; Fume Hood; PCR Cabinet; Dispensing Booth.

Date of issue: July 22, 2022

Date of expiry: July 21, 2025

Date of change: October 20, 2022

General Manager:

**BEIJING HUAGUANG CERTIFICATION
OF MEDICAL DEVICES CO., LTD.**



BIOBASE®

Add: No.9 Gangxing Road, High-tech Zone, Jinan City,
shandong Province, China

Tel: +86-531-81219803/01 Fax: +86-531-81219804

Email: export@biobase.com

Website: www.biobase.cc www.meihuatrade.com



BIOBASE CHINA



BIOBASE CHINA



BIOBASE CHINA