



# SKH029

Carro de cirugía





## Control de calidad importante

- \* Cuerpo completo de acero inoxidable.
- \* Dos estantes de acero inoxidable, dos cajones con asa de plástico y con cerradura, tamaño del cajón (LxAxAI): 100x280x350mm
- \* Empuje la manija en el lado derecho
- \* Ruedas cubiertas de 3", dos con freno
- \* Barandillas en tres lados en cada capa
- \* Marco de tubo redondo,  $\Phi 25 \times 1.2 \text{mm}$  y  $\Phi 13 \times 1 \text{mm}$

## Parámetros técnicos

Tamaño L*A*A	670*450*900mm
--------------	---------------

## Configuración técnica

Tary y cubo de basura	1 pieza
Cerradura de 25mm para carro	2 piezas
Carril de cajón de 24"	1 par



## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N 2 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
international City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue