




X12


Mesa de examinación
eléctrica

SKD-A



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098



Altura ajustable

Parámetros técnicos

External L*W*H 1900*600* (500-800) mm

Configuración técnica

Motor lineal	1 pc
Cubrepiés antideslizante	4 pcs
Almohada extraíble	1 pc

Electrónica

Voltaje	100-240V
Frecuencia	50/60Hz
Enchufe	as required (hospital grade)

QUALITY CONTROL

Cuero

Ignífero, resistente al amarilleamiento, resistente al moho, resistente a la flexión 100.000 veces resistente al desgaste

100,000

España

Ignífero, 4cm de espesor, espuma de alta densidad, densidad 30kg/m

4 | 30

Pintura

ASTM pruebas anti-material
Espesor de la pintura 0.12mm, brillo 60, la pintura puede resistir un impacto de 50kg

0.12 | 60 | 50

Material metálico

Baosteel en la lista Fortune Global 500
Panasonic robótica

500 | 360°

Sistema eléctrico

Taiwán motor
Inferior a 50dB, con vida útil 20.000 veces de movimiento

50 | 20,000

Piezas ricas



1.2-1.5 mm Espesor del bastidor de acero con recubrimiento de polvo, 250kg de carga portante



4 fundas antideslizantes para pies de cama



Controlador de mano, mantener calificado después de 2000 veces de uso





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC.
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.