

Removable Double Nozzle Eye Wash

EWM-201



Specifications

Model	EWM-201
Body	SS304
Function	Eye and face washing
Surface finishing	Sanbed
Spray head feature	Patented design, each spray head equips with one water-flow control valve
Dust cover	Integrated with the spray head (not connected by chain)
Logo	The logo of the switch is placed on the front of the push board,which indicate the direction of switch, which is eye-catching and clear.
Push board	Push the board by hand, you can wash your eyes and face, easy to use
ball valve	SS304
Working pressure	0.2-0.4Mpa
Installation	Desktop

Bench Eye Wash

EEW-100



Description

- With cold and hot water mixing function, it has accident flushing and normal water supply in the laboratory.
- The washing part is equipped with water supply hose to extend the washing range.
- The water supply elbow can be rotated 360 times for easy operation.

Specifications

Model	EEW-100
Body	Solid brass
Function	eye and face washing
Coating	High gloss thick plating coating, resistant to chemicals,UV fading and heat
Spray head	After slow pressure treatment, the water will be foamy to avoid secondary damage
Dust cap	PP materials,automatically turned on by water when used
Push board	push the board by hand, you can wash your eyes and face, easy to use
ball valve	SS304
Working pressure	0.2-0.4Mpa