



SKP011

SECTIONAL MEMORY
FOAM MATTRESS



SAIKANG / Medical

Healthier
your life

SKP011

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 1925×900×140mm

Zippered Cover

- Material PU
- Color grey or dark blue
- Size (LxWxH) 1925×900×140mm

Upper Layer

- Material visco-elastic memory foam
- Color blue, green, black
- Thickness 40mm
- Density 42-45kg/m³

Firm Base Layer

- Material polyurethane foam
- Color creamy-white
- Thickness 100mm
- Density 40kg/m³

Upper Layer

- a) High density: 42-45kg/m³
- b) Unique thermo active properties to provide good support and excellent pressure relief



Firm Base Layer

- a) High density: 40kg/m³
- b) Very elastic and thermosetting. Sectional design works to make the mattress and bed surface fit well in different positions such as back-rest and knee-rest adjustment.



Zippered Cover

- Air through, water proof, soft, wind-resistant, strong tensile strength and anti-infection. Long zipper design, easy to remove and wash.



Size

- Customized size acceptable





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue