



SKP012

18 Cells Tubular Air Mattress



SKP012

18 Cells Tubular Air Mattress

Pump

- Power supply: 110V/220V/50-60Hz
- Air output: 6-8 liter/min.
- Pressure range: 50mmHg - 120 mmHg
- Cycle time: 12 min.
- Power consumption: 7 W
- Fuse: 1A
- Size: 27.5 (L) x 13.5 (W) x 9 (H) cm
- Weight: 1.6kg

Mattress

- Standard Size: 80 "(L) x 32-36 "(W)
- 4" 24 cells , 4.5" 18 cells
- Top cover: Nylon PU /strength
- Air cells: Nylon / PVC/TPU
- Bottom: tarpaulin /nylon PU
- Mattress weight : 5.5-6.5kg
- Weight Capacity: 135-145kg