



# SKB039(E)

## AMBULANCE STRETCHER TROLLEY

Registration code: SKC-3F5



**SAIKANG** / Medical

Healthier  
your life

# SKB039(E)

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ Full size (LxWxH)	1900x540x800mm
▪ Lowest size (LxWxH)	1900x540x200mm
▪ Back-rest lifting	0-85°
▪ Net weight	45kg
▪ Load bearing	≤159kg

#### Side rail

Side rail is foldable.



#### Safety lock

With safety lock, it must be locked up when transport the patients, and opened out when get on the ambulance.



#### Caster

5' silence caster. Advanced protection universal brake wheel, easy to move, simple to operate, reliable performance.



#### Function

Can be folded into a chair, the angle of the chair can be adjusted in five gear.



#### Fixing devices

Can be locked by the fixing devices when on the ambulance car.



#### Structure

Premium alluminium alloy structure, more stable, but light weight.



#### Leather / Safety belts

Artificial leather, bright appearance, soft texture, abrasion resistance, folding resistance, acid and alkali resistance. With 2 safety belts, make patient safer.



#### Sponge

Flame retardant, 5 cm thickness, high-density foam, density 30 kg/m<sup>3</sup>.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue