

Water Bath

WB-1R1H-3 WB-1R2H-7 WB-2R4H-15 WB-2R6H-25



Features

AI Intelligent temperature controller with PID program, digital tube display and function of timing and over-temperature protection.

U Add the electric drain switch, one-button operation, convenient to drain, Note: single hole water baths have no drain function.

Shield Power off when water shortage to ensure the safe operation of the product.



Stainless steel stretching inner container and roof cover.



Specifications

Model	WB-1R1H-3	WB-1R2H-7	WB-2R4H-15	WB-2R6H-25	
Classification Specifications	1 hole	1 row & 2 holes	2 rows & 4 holes	2 rows & 6 holes	
Heating Mode	Nature water convection heat transfer				
Function	Temp. Range				
	Temp. Resolution Ratio				
	Temp. Motion				
	Temp. Uniformity				
Structure	Inner Chamber				
	Outer Shell				
	Heater				
	Power Rating				
Controller	Temp. Control Mode				
	Temp. Setting Mode				
	Temp. Display Mode				
	Timer				
	Operation Function				
	Program Mode				
	Additional Function				
Sensor					
Safety Device					
Specification	Inner Chamber size(W*L*H)(mm)	150*135*150	300*135*150	325*300*150	500*300*150
	Exterior size (W*L*H)(mm)	170*154*210	318*168*210	350*318*210	524*322*210
	Packing size (W*L*H)(mm)	240*230*280	390*240*280	420*390*280	600*390*280
	Volume	3.4L	6.1L	14.6L	22.5L
	Load per rack	5kg			
	Shelf Number	1			
	Power Supply (50/60Hz)	AC220V/1.8A	AC220V/2.3A	AC220V/4.5A	AC220V/6.8A
Current rating NW/GW (kg)	3.3/4	4.5/5	6/7	7.5/9	