

# Infitek

---



# INTRODUCTION

Founded in 2010, with the mission of continuously improving the intelligence, precision, safety and convenience of the laboratory, Infitek is a professional manufacturer driven by innovation and service in laboratory and Medical field, our company is certified by ISO9001, ISO13485 and Intellectual Property Management System

In addition to providing our self-produced products, we also provide customers with purchasing and inspection services in whole China market. We provide fully after-sales service to the products running through us.



Each of our facilities have more than 10 years of experienced R&D personnel to ensure the continuous improvement of product quality. In order to ensure processing accuracy, we use internationally renowned sheet metal processing equipment and CNC machine tools, such as DMG MORI and AMADA. The whole production process is strictly implemented in accordance with industry standards. All accessories, metal processing, assembly and inspection work are completed by highly trained and responsible employees.


To guarantee the performance of the finished product, each product line has set up an independent testing room.

**18** Factories

We have 18 factories in form of self-establishment, holding or equity participation.



# CERTIFICATE



ISO 13485




ISO 9001



FDA



Intellectual Property Management System Certificate



Enterprise Credit Evaluation : AAA

# Stainless Steel Water Distiller Water Control Type



WDST-5    WDST-10    WDST-20

## Features

- The whole machine is made of high-quality stainless steel.
- Coil type condensing structure, good cooling effect, large water output.
- The cooling water flows back into the evaporation barrel to save water resources.
- The drain valve can drain the concentrated water in the barrel to improve fouling conditions.
- The heating power is automatically cut off when the water is lacking, and the heating will continue after the return water is replenished to the working water level.
- The copper immersion heating tube has high thermal efficiency and long service life.



## Specifications

Model	WDST-5	WDST-10	WDST-20
Applicable power supply	220V/50/60Hz	380V/50Hz	380V/50Hz
Power	4.5KW	7.5KW	15kW
Water output	≥5L/h	≥10L/h	≥20L/h
Packing size (wood) (cm)	35×34×78	42×38×90	51×36×80



Infitek



Infitek



Infitek



Infitek



Infitek

# Infitek

## Infitek Co., Ltd.

TEL: +86-531-88982330

FAX: +86-531-88983691

Website: [infitek.com](http://infitek.com)

Email: [info@infitek.com](mailto:info@infitek.com)

Service: [support@infitek.com](mailto:support@infitek.com)

Address: Rm. 2014, Bldg. 3, Ligaoguoji Huayuan, No. 1222, West Aoti Road, Lixia District, Jinan, Shandong

## US Office

### INFITEK INC.

522W RIVERSIDE AVE STE N, SPOKANE, WA 99201

EMAIL: [INFO@INFITEK.COM](mailto:INFO@INFITEK.COM)