

Table Top Autoclave

Class N Series



Introduction:

Class N autoclave discharges the cold air in the sterilization chamber by gravity displacement, saturates the internal steam, and reaches the set pressure and temperature. To achieve the purpose of sterilizing the items to be sterilized. Suitable for stomatology, ophthalmology, beauty, operating room and other places. It is suitable for sterilization of glassware, bare metal utensils, and unencapsulated solid instruments.

Features:

- *High quality stainless steel sterilization chamber and sealing cover.
- *Microcomputer control for automatic circulation, with 121°C and 134°C freely settable.
- *LED display with touch buttons.
- *Buzzer reminder at the end of sterilization
- *Automatically discharge cold air. Automatically discharge steam at the end of sterilization.
- *Built-in heating tube device to quickly provide saturated steam.
- *Two sterilization procedures at 121°C and 134°C are selectable.

Technical Parameters:

Model	BKX-T18E	BKX-T23E
Capacity	18L	23L
Door	Manual	
Display	LED	
Max Designed Pressure	0.25MPa	
Working Pressure	0.22MPa	
Max Design Temp.	139°C	
Working Temp.	121/134°C	
Temp. Precision	0.1°C	
Chamber Material	SUS304	
Chamber Size	Φ 250*370mm	Φ 250*452mm
Altitude	≤2500m	
Rack	1 rack with 2 trays	
Standard Accessories	Protective gloves, draining hose, door seal	
Consumption	1.5KW	2KW
Power Supply	AC220V, 50Hz (Standard) ; AC110V, 60Hz (Optional)	
External Size(W*D*H)	500*400*350mm	600*400*350mm
Net Weight	24kg	27kg
Package Size(W*D*H)	540*450*390mm	640*450*390mm
Gross Weight	26kg	29kg